











More Than Just A Dental Lab. . .

# We Decided To Make It Clearer

This overview of available DTS Crown and Bridge restorations lets you compare all the modern and traditional restorations we make. On this page we have carefully chosen the restorations we currently recommend and believe are all you need to complete quality dentistry everyday.

	 <b>VITA SUPRINITY®</b>	 <b>Core Zirconia™</b>	 <b>Opalite™</b>	 <b>Essential™ Range</b>
<b>Material</b>	Zirconia reinforced Lithium Silicate Glass Ceramic	Zirconia coping with porcelain overlay	Monolithic translucent zirconia	Non-precious bonded
<b>Indications</b>	Anterior Crowns, Posterior Crowns	Single unit crowns, multi unit bridges, full arch	Bruxers/grinders, implants	Single unit crowns, multi unit bridges
<b>Recommended Use</b>	Crowns, inlays/onlays, veneers, crowns on implant abutments	Full mouth Crown and Bridge	Mostly for Posteriors, where high strength required	Full mouth
<b>Preparation Guide</b>	 Recommended Preparation For Anterior Teeth    Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth    Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth    Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth    Recommended Preparation For Posterior Teeth
<b>DTS Recommends</b>	Recommended for single anterior or posterior restorations where the highest level of aesthetics are to be achieved	Recommended for multiple crown and bridge cases where high aesthetics are to be achieved	Recommended for posterior crown and bridge cases only where strength is of the utmost importance, ideal for natural looking replacement to metal full crowns	Recommended for anywhere in the mouth in crown and bridge cases where cost is the highest priority
<b>Cementation</b>	Conventional cementation or adhesively bonded	Conventional cementation or adhesively bonded	Conventional cementation or adhesively bonded	Conventional cementation
<b>Turnaround Time</b>	5 Lab Days	5 Lab Days	5 Lab Days	10 Lab Days
<b>Aesthetics</b>	✓✓✓✓✓	✓✓✓✓✓	✓✓✓	✓✓✓
<b>Strength</b>	✓✓✓✓	✓✓✓✓	✓✓✓✓✓	✓✓✓✓
<b>Warranty</b>	5 years	5 years	5 years	1 year

dts-international.com

					
	Standard PFM	Lava™ Classic	IPS e.max® CAD	Lava™ Ultimate	VITA ENAMIC
Material	Semi-precious or precious alloy coping with layered porcelain overlay	Independent 3M zirconia coping, with multi-layered, porcelain overlay	Ivoclar e.max lithium disilicate	Nano-composite	Ceramic strengthened by a polymer network
Indications	Single unit crowns, multi-unit bridges	Full Mouth, implants	Anterior and posterior crown	Single unit, implants	Anterior Crowns, Posterior Crowns
Recommended Use	Full mouth	Full mouth Crown and Bridge	Full Mouth Crowns, Onlays/ Inlays, veneers	Crowns, inlays/onlays, veneers, recommended posterior implant	Crowns, inlays/onlays, veneers, recommended posterior implant
Preparation Guide	 Recommended Preparation For Anterior Teeth      Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth      Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth      Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth      Recommended Preparation For Posterior Teeth	 Recommended Preparation For Anterior Teeth      Recommended Preparation For Posterior Teeth
Cementation	Conventional cementation	Conventional cementation or adhesively bonded	Conventional cementation or adhesively bonded	Adhesive bonding	Adhesive bonding
Turnaround Time	5 Lab Days	5 Lab Days	5 Lab Days	5 Lab Days	5 Lab Days
Aesthetics	✓✓✓	✓✓✓✓✓	✓✓✓✓✓	✓✓✓	✓✓✓
Strength	✓✓✓✓	✓✓✓✓	✓✓✓	✓✓✓	✓✓✓✓
Warranty	1 year	5 years	5 years	5 years	5 years

For Further Information On Our Full Range Of Products Or Services:

✉ 8 Annfield Place | Duke Street | Glasgow | G31 2XQ

☎ Customer Services: +44 141 556 5619 | Collection Line: +44 141 554 5657

@ restorations@dts-international.com



To find out more about Digital Dentistry and how it can benefit your practice, email [digital@dts-international.com](mailto:digital@dts-international.com)